







HME Cap™

Disclaimer

Atos Medical offers no warranty - neither expressed nor implied - to the purchaser hereunder as to the lifetime of the product delivered, which may vary with individual use and biological conditions. Furthermore, Atos Medical offers no warranty of merchantability or fitness of the product for any particular purpose.

Patents and trademarks

Provox® is a registered trademark owned by Atos Medical AB, Sweden. Provox® HME Cap™ is a trademark of Atos Medical AB.

For information about protective rights (e.g. patents), please refer to our web page www.atosmedical.com/patents.

Manufacturer; Hersteller; Fabrikant; Fabricant; Produttore; Fabricante; Fabricante; Tillverkare; Producēt; Produsent; Valmistaja; Framleiðandi; Tootja; Ražotājs; Gamintojas; Výrobce; Gyártó; Výrobca; Proizvajalec; Producent; Producător; Proizvodčač; Proizvodčač; Катаскенасстής; Производител; Üretici; მწარმოებელი; Արտադրող; İstehsalçı; Производитель; Produsen; Pengilang; 製造業者; 제조사; 製造商; 制造商; היצרן; المصنعة; الجهة المصنعة

Date of manufacture; Herstellungsdatum; Productiedatum; Date de fabrication; Data di fabbricazione; Fecha de fabricación; Data de fabrico; Tillverkningsdatum; Fremstillingsdato; Produksjonsdato; Valmistuspäivä; Framleiðsludagsetning; Valmistamise kuupäev; Ražošanas datums; Pagaminimo data; Datum výroby; Gyártás időpontja; Dátum výroby; Datum proizvodnje; Data produkcji; Data fabricației; Datum proizvodnje; Datum proizvodnje; Ημερομηνία κατασκευής; Дата на производство; Üretim tarihi; წარმოების თარიღი; Արտադրման ամսաթիվ; İstehsal tarixi; Дата изготовления; Tanggal produksi; Tarikh pembuatan; 製造日; 제조일자; 製造日期; 制造日期; تاريخ التصنيع; תאריך ייצור

Single Patient – multiple use; Einzelner Patient – mehrfach anwendbar; Eén patiënt – meervoudig gebruik; Un seul patient – à usage multiple; Singolo paziente – uso multiplo; Un solo paciente – uso múltiple; Paciente único – várias utilizações; En patient – flera användningar; Enkelt patient – flergangsbrug; Kun til bruk på én pasient – flergangsbruk; Potilaskohtainen – voidaan käyttää useita kertoja; Fyrir einn sjúkling – margnota; Ühel patsiendil korduvalt kasutata; Vi en spaci en ts – vairākkārtēja lietošana; Vienas pacientas – daugkartinis naudojimas; Jeden pacient – vícenásobné použití; Egyetlen beteg esetében többször újra használható; Jeden pacient – viacnásobné použitie; En bolnik – večkratna uporaba; Wielokrotne użycie u jednego pacjenta”; Un singur pacient – utilizare multiplă; Jedan pacijent – višestruka uporaba; Један пациент –

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Instructions for use

Intended use

Device for rehabilitation after total laryngectomy.

The Provox HME Cap is a dome-shape titanium ring that allows you to use a special Provox FreeHands HME cassette (REF 8220, 8221) without a FreeHands valve unit. It cannot be used with any other type of HME cassette.

Attach the device by means of any Provox Adhesive base plate or Provox LaryTube. Close the front opening of the cap manually when speaking with a voice prosthesis.

Replace the HME cassette every 24 hours (more frequently if needed). The cap itself can be cleaned and reused.

Because the device is located in front of your tracheostoma and thus your airways, please observe the following warnings:

Warnings

- **ALWAYS assemble the device (HME cassette and cap) carefully before use.**
- **If the cap should detach from the HME cassette during use, you must remove the complete adhesive base plate from the stoma, with the HME cassette still mounted in the base plate. Do not try to manipulate in any other way, and NEVER try to re-attach the cap.**
- **NEVER remove the foam from the HME cassette. This would eliminate the HME features and can damage the cassette casing. It can also increase the risk of inhaling particles from the device.**

- **Do NOT reuse the HME cassette or the foam, e.g. after washing out. Washing eliminates the calcium chloride needed for a sufficient HME function and can increase the risk of infection.**

CAUTION: Do not use the device until it is completely dry. Inhalation of disinfectant fumes can cause severe coughing and irritation of the airways.

If you have further questions, please contact our customer service department.

Disposal

Always follow medical practice and national requirements regarding biohazard when disposing of a used medical device.

Reporting

Please note that any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the national authority of the country in which the user and/or patient resides.

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