TRACOE CONNECTION TUBE WITH CHECK VALVE

REF ODG0702SC167

EN	INSTRUCTION FOR USE
DA	BRUGSANVISNING
DE	GEBRAUCHSANWEISUNG
ES	INSTRUCCIONES DE USO
FR	MODE D'EMPLOI
IT	ISTRUZIONI PER L'USO
NO	BRUKSANVISNING
SV	BRUKSANVISNING

 Promepia S.A.M.
9 Avenue Albert II - 98000 MONACO

 Coloplast A/S
Høtvedam 1
3050 Humlebaek
DENMARK

 Promepia Group S.à.r.l.44 Avenue J-F Kennedy 1855
Luxembourg

REF: N039-0911
REV: 10
DATE: 09/2025

EN

Instruction for Use

TRACOE Connection Tube with check valve (ODG0702SC167) for connecting the TRACOE Cuff Pressure Monitor Sensitive and TRACOE Smart Cuff Manager to a tracheostomy tube or an endotracheal tube.

1. **Pack contents**
The package of the TRACOE Connection Tube with check valve (ODG0702SC167) contains 10 units individually packed connection tubes for the TRACOE Cuff Pressure Monitor Sensitive or TRACOE Smart Cuff Manager.

All tubes are fitted with a female Luer connector with check valve at one end and a male Luer lock connector at the other end.
Length of the connection tubes: 1 meter

2. **Intended Use**
The TRACOE Connection Tube is designed for the pneumatic connection between the TRACOE Cuff Pressure Monitor Sensitive/TRACOE Smart Cuff Manager and the valve of the cuff inflation line/pilot balloon of a tracheostomy tube or an endotracheal tube.
Duration of use: The maximum duration of use is 72 hours in total. Products that show only minor damage must be replaced.

3. **Indication**
This device is intended for use in adults, pediatrics, and neonates intubated by tracheostomy tube or an endotracheal tube. Multiple use on the same patient is allowed. This device is intended for use by qualified and trained personnel.

4. **Instructions**
To create a pneumatic connection between the TRACOE Cuff Pressure Monitor Sensitive or the TRACOE Smart Cuff Manager and the cuff of a tracheostomy tube or an endotracheal tube, the TRACOE Connection Tube with its female Luer connector with check valve is pushed firmly onto the corresponding male component on the TRACOE Cuff Pressure Monitor Sensitive or the TRACOE Smart Cuff Manager. The male Luer connector located on the other end of the tube is pushed or screwed firmly onto the valve of the cuff inflation line/pilot balloon supplying the cuff on the tracheostomy tube or endotracheal tube.
After establishing the connection and inflating the cuff, please monitor the pressure stability for a minute using the TRACOE Cuff Pressure Monitor or TRACOE Smart Cuff Manager. A stable pressure maintenance indicates a secure seal between the cuff inflation line/pilot balloon of the tracheostomy tube or an endotracheal tube and the TRACOE Connection Tube.

5. **Contraindications**
Not for use as a standalone airway device.

6. **Cautions/Warnings**
• Do not reuse on another patient.
• Do not use if the packaging is damaged.
• Do not use if it has exceeded the expiry date.
• Check the integrity of the devices before any use. Do not use if there are any suspicions regarding integrity.
• Ensure compatibility with the tracheostomy tube and respiratory circuit.
• Regularly monitor the connection for disconnection or secretion buildup.
• Use only under supervision of qualified healthcare personnel.
• Not intended for sterilization.
• Ensure secure fixation to avoid accidental disconnection.

7. **Reporting serious incidents**
Any serious incident related to the use of the TRACOE Connection Tube must be immediately reported to the manufacturer and the competent authority of the Member State in which the user or patient is based. These incidents may include (but are not limited to) serious injuries, worsening of state of health or potential safety risks.

8. **Conditions of storage & Use**
Storage and transport in a normal environment, i.e. protected from extreme temperatures and humidity.

9. **Disposal**
After use, get rid of the TRACOE Connection Tube and its packaging in accordance with the current protocols in the health establishment, as well as with the local and national regulations, and/or instructions of the competent authorities.

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