

Justification:	Function:	Electronic signature justification:	Signed: Date (yyyy-mm-dd) - Time (hh:mm):
Issued:	QA	Carolina Johansson - SEHRBJNC	2021-12-06 - 09:59
Reviewed:	QA	Karolina Nilsson - KARNIL	2021-12-07 - 09:56
Approved:	DD	Diana Tieger - DIATIE	2021-12-07 - 10:35
Released:	QA	Carolina Johansson - SEHRBJNC	2021-12-27 - 14:37

This document has been electronically signed by the persons above.

# Product Information

## Provox Cleaning Towel



**Product description:**

Provox Cleaning Towel is alcohol-based and non-perfumed and is available in resealable pouches.

Document No: 10000042798 Edition: 02 Release date: 2021-12-27

Released

## Product Information

<b>Document ID:</b>	PF076-03-TechInfo	<b>Edition:</b>	02
<b>Manufacturer:</b>	Atos Medical AB Kraftgatan 8 SE-242 35 Hörby, Sweden		
<b>Classification: (EU) 2017/745</b>	Class I, Rule 1		
<b>Intended Use:</b>	Provox Cleaning Towel is intended for cleaning around the stoma, it will remove oil from the skin. They are indented to use before application of Provax Adhesives.		
<b>Use specifications:</b>	N/A		
<b>Contraindications:</b>	N/A		
<b>CE Mark:</b>	Yes. Devices are CE-marked.		
<b>GMDN code:</b>	46205		
<b>Sterilization:</b>	Non-sterile		
<b>Raw material:</b>	Towel: Spunlaced viscose, Polyester Solution: Chlorhexidine, Ethanol, BTC 2125, deionised water. MSDS-7244-SE		
<b>Latex information:</b>	Not manufactured with natural rubber latex		
<b>Biological origin:</b>	The device is not manufactured with materials derived from human or animal source.		
<b>Handling and storage:</b>	Store the product dry and away from sunlight at room temperature. Excursions permitted between 2°C - 42°C.		
<b>Waste handling and disposal:</b>	Waste handling and disposal should be carried out in agreement with medical practice and applicable national laws and legislations. Used product may be a potential biohazard.		
<b>Hazardous components:</b>	The product contains ethanol (concentrated ethanol is flammable) and chlorhexidine (concentrated chlorhexidine may cause eye damages if in contact with the eyes).		
<b>Expiration date:</b>	2 years after manufacturing.		
<b>Packaging:</b>	10 towels are packed together in a plastic bag made of polyester/polyethene and PT 12 (white)/PE 70. 20 plastic bags are packed together in every cardboard box		

## Product Information

Devices under Basic UDI-DI: 7331791-ADH-A-000-0003-UH

REF	Name	UDI-DI
7244*	Provox Cleaning Towel 10-p	07331791001536
7244**	Provox Cleaning Towel 10-p	07331791015762
7244-18	Provox Cleaning Towel	07331791014758

\*Contains 20 bags with 10 towels

\*\*Contains 1 bag with 10 towels

### Atos Medical AB compatible products:

Range	BASIC UDI-DI
Provox Luna Adhesive	7331791-ADH-0-000-0000-CQ
Provox StabiliBase	7331791-ADH-0-000-0000-CQ
Provox Adhesive Flexiderm	7331791-ADH-0-000-0000-CQ
Provox Adhesive Optiderm	7331791-ADH-0-000-0000-CQ
Provox XtraBase Adhesive	7331791-ADH-0-000-0000-CQ
Provox Skin Barrier	7331791-ADH-A-000-0004-UL
Provox Adhesive Remover	7331791-ADH-A-000-0005-UP
Provox Silicone Glue	7331791-GEN-A-000-0003-EF
Provox Life Standard Adhesives	7331791-ADH-0-000-0001-CT
Provox Life Sensitive Adhesives	7331791-ADH-0-000-0001-CT
Provox Life Stability Adhesive	7331791-ADH-0-000-0001-CT
Provox Life Night Adhesive	7331791-ADH-0-000-0001-CT